

ITGAV MAB (ABITUZUMAB)

货号: 28138

产品全名: ITGAV(Abituzumab) 单抗

基因符号 Integrin alpha-V;CD51

描述: ITGAV(abituzumab) 单抗

背景: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

经过测试的应用: Flow Cyt

推荐稀释比: Flow Cyt 1:100

种属反应性: Humanized

亚型: IgG2

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: 人 ITGAV

成分: Lyophilized from sterile PBS, pH 7.4. 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form.

After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. ITGAV protein is highly expressed on the surface of A549 cell membrane. Flow cytometry analysis with $1\ \mu\text{g}/\text{mL}$ Anti-ITGAV(abituzumab biosimilar) mAb (28138) (Blue histogram) or isotype control mAb (Red histogram) on A549 cells.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
