

## HUMAN ACKR3 FULL LENGTH PROTEIN

货号: 11076

产品全名: 人 ACKR3 Full Length 蛋白

规格: 10/50/100 µg

基因符号 CMKOR1; CXC-R7; CXCR-7; CXCR7; GPR159; RDC-1; RDC1

目标蛋白: ACKR3

**UNIPROT ID:** P25106

描述: 人 ACKR3 full length protein-synthetic nanodisc

背景: This gene encodes a member of the G-protein coupled receptor family. Although this protein was earlier thought to be a receptor for vasoactive intestinal peptide (VIP), it is now considered to be an orphan receptor, in that its endogenous ligand has not been identified. The protein is also a coreceptor for human immunodeficiency viruses (HIV).

Translocations involving this gene and HMGA2 on chromosome 12 have been observed in lipomas. [provided by RefSeq, Jul 2008]

物种/宿主: HEK293

分子量: The human full length ACKR3 protein has a MW of 41.3 kDa

**Formulation & Reconstitution:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form.

After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.