

## ANGPTL3 (DM205) RABBIT MAB

货号: 28548

产品全名: ANGPTL3(DM205) 兔单克隆抗体

基因符号 ANG-5; ANGPT5; ANL3; FHBL2

描述: ANGPTL3 antibody(DM205) 兔单克隆抗体

背景: This gene encodes a member of a family of secreted proteins that function in angiogenesis. The encoded protein; which is expressed predominantly in the liver; is further processed into an N-terminal coiled-coil domain-containing chain and a C-terminal fibrinogen chain. The N-terminal chain is important for lipid metabolism; while the C-terminal chain may be involved in angiogenesis. Mutations in this gene cause familial hypobetalipoproteinemia type 2.

经过测试的应用: ELISA; Flow Cyt

推荐稀释比: ELISA 1:5000-10000; Flow Cyt 1:100

种属反应性: Rabbit

亚型: Rabbit IgG

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: 人 ANGPTL3

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form.

After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with Anti-ANGPTL3 (DM205) on Expi293 cells transfected with human ANGPTL3 (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).