

15 PGDH RABBIT MAB

货号: N261779

产品全名: 15 PGDH 兔单克隆抗体

基因符号 15-PGDH; Hpgd; PGDH; PGDH1; PHOAR1; SDR36C1

UNIPROT ID: P15428

背景: Prostaglandin inactivation. Contributes to the regulation of events that are under the control of prostaglandin levels. Catalyzes the NAD-dependent dehydrogenation of lipoxin A4 to form 15-oxo-lipoxin A4. Inhibits in vivo proliferation of colon cancer cells.

抗原: Recombinant protein of human 15-PGDH

经过测试的应用: WB,IHC-F,IHC-P,ICC/IF

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R09-3D2

分子量: Calculated MW: 29 kDa; Observed MW: 25 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Immunology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of 15PGDH in SW480 lysates using 15 PGDH antibody.