

## ABAT RABBIT MAB

货号: N261783

产品全名: ABAT 兔单克隆抗体

基因符号 ABAT; GABA transaminase; GABA transferase; GABAT; LAIBAT

**UNIPROT ID:** P80404

背景: Catalyzes the conversion of gamma-aminobutyrate and L-beta-aminoisobutyrate to succinate semialdehyde and methylmalonate semialdehyde, respectively. Can also convert delta-aminovalerate and beta-alanine.

抗原: A synthetic peptide of human ABAT

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R08-9E3

分子量: Calculated MW: 56 kDa; Observed MW: 56 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Rat

成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Neuroscience

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of ABAT/GABAT in rat Brain lysates using ABAT antibody.

