

ACADS RABBIT MAB

货号: N261785

产品全名: ACADS 兔单克隆抗体

基因符号 SCAD; ACAD3

UNIPROT ID: P16219

背景: Introduces a double bond at position 2 in saturated acyl-CoA's of short chain length, i.e. less than 6 carbon atoms.

抗原: A synthetic peptide of human ACADS

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R08-4C2

分子量: Calculated MW: 44 kDa; Observed MW: 44 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of ACADS in Jurkat, C6, HeLa lysates using ACADS antibody.



Western blot analysis of ACADS in mouse liver, 3T3 lysates using ACADS antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
