

ACETYL-HISTONE H4 (LYS16) RABBIT PAB

货号: N225522

产品全名: Acetyl-Histone H4 (Lys16) 兔多抗

基因符号 H4K16ac; Histone 1 H4 antibody; HIST1H4K antibody; MGC24116 antibody; dJ160A22.1 antibody

UNIPROT ID: P62805

背景: The chromatin fiber is further compacted through the interaction of a linker histone, H1, with the DNA between the nucleosomes to form higher order chromatin structures. This gene is intronless and encodes a member of the histone H4 family.

抗原: Synthetic Peptide of Histone H4 (Acetyl Lys16)

经过测试的应用: ICC/IF, WB, IHC-F, IHC-P

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 11 kDa; Observed MW: 14 kDa

亚型: IgG

纯化: Affinity Purified


种属反应性: Human, Mouse and Rat

Modification: Acetylated

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics and Nuclear Signaling

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Acetyl-Histone H4 (Lys16) in HeLa lysates untreated or treated with TSA lysates using Acetyl-Histone H4(Lys16) antibody(upper) and betaactin (5B7-) antibody(lower) TDY041.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
