

ACETYL-LYSINE RABBIT PAB

货号: N226080

产品全名: Acetyl-Lysine 兔多抗

基因符号 -

UNIPROT ID: -

背景: Acetylation of lysine residues within proteins has emerged as an important mechanism used by cells to overcome this repression.

抗原: Acetylated BSA

经过测试的应用: WB

推荐稀释比: WB: 1/500-1/1000

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Refer to figures

亚型: IgG

纯化: Affinity Purified


种属反应性: Species independent

Modification: Acetylated

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics and Nuclear Signaling

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Acetyl-Lysine in HeLa lysates using Pan Acetyl-Lysine antibody.