

ACETYL-P53 (LYS382) (2G4) MOUSE MAB

货号: N261359

产品全名: Acetyl-p53 (Lys382) (2G4) 小鼠单抗

基因符号 Cellular tumor antigen p53; Cys 51 stop; Tp53; Tumor protein p53

UNIPROT ID: P04637

背景: Tumor protein p53, a nuclear protein, plays an essential role in the regulation of cell cycle, specifically in the transition from G0 to G1. It is found in very low levels in normal cells, however, in a variety of transformed cell lines, it is expressed in high amounts, and believed to contribute to transformation and malignancy. p53 is a DNA-binding protein containing DNA-binding, oligomerization and transcription activation domains.

抗原: Synthetic peptide conjugated to KLH.

经过测试的应用: IHC-P

推荐稀释比: IHC: 1/50-1/100

种属反应性: Mouse

克隆性: Mouse Monoclonal

克隆编号: 2G4-5C4-4D2

分子量: -

亚型: IgG1

纯化: Affinity Purified


种属反应性: Human,Rat,Mouse

Modification: Acetylated


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemical analysis of paraffin-embedded Human tonsils using Acetyl-p53 (Lys382) (2G4) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human Breast Carcinoma Tissue using Acetyl-p53 (Lys382) (2G4) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.