

## ACETYL-SMAD2/3 (LYS19) RABBIT PAB

货号: N225507

产品全名: Acetyl-Smad2/3 (Lys19) 兔多抗

基因符号 JV15-2; SMAD 3; hMAD-2; SMAD 2

**UNIPROT ID:** Q15796

**背景:** The protein encoded by this gene belongs to the SMAD, a family of proteins similar to the gene products of the Drosophila gene 'mothers against decapentaplegic' (Mad) and the C. elegans gene Sma. SMAD proteins are signal transducers and transcriptional modulators that mediate multiple signaling pathways.

**抗原:** Synthesized peptide derived from human Smad2/3 (Acetyl-Lys19)

**经过测试的应用:** IHC-P,WB

**推荐稀释比:** WB: 1/500-1/1000 IHC: 1/50-1/100

**种属反应性:** Rabbit

**克隆性:** Rabbit Polyclonal

**分子量:** Calculated MW: 52 kDa; Observed MW: 60 kDa

**亚型:** IgG

**纯化:** Affinity Chromatography

**种属反应性:** Human, Mouse and Rat

**Modification:** Acetylated

**成分:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**研究领域:** Signal Transduction

**储存和运输:** Store at -20°C. Avoid repeated freezing and thawing