

ACTIVIN RECEPTOR TYPE IA RABBIT MAB

货号: N261795

产品全名: Activin Receptor Type IA 兔单克隆抗体

基因符号 ACTRI; Acvr1; ACVRIA; ACVRLK2; ALK2; FOP; SKR1; TSRI

UNIPROT ID: Q04771

背景: On ligand binding, forms a receptor complex consisting of two type II and two type I transmembrane serine/threonine kinases. Type II receptors phosphorylate and activate type I receptors which autophosphorylate, then bind and activate SMAD transcriptional regulators.

抗原: A synthetic peptide of human Activin Receptor Type IA

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R07-7G4

分子量: Calculated MW: 57 kDa; Observed MW: 57 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human,Mouse


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Activin Receptor Type IA in Hela lysates using Activin Receptor Type IA antibody.



Western blot analysis of Activin Receptor Type IA in mouse heart, mouse brain lysates using Activin Receptor Type IA antibody.