

## ADSS 2 RABBIT MAB

货号: N261805

产品全名: AdSS 2 兔单克隆抗体

基因符号 ADEH; ADSS 2

**UNIPROT ID:** P30520

**背景:** Plays an important role in the de novo pathway and in the salvage pathway of purine nucleotide biosynthesis. Catalyzes the first committed step in the biosynthesis of AMP from IMP.

**抗原:** Recombinant protein of human AdSS 2

**经过测试的应用:** WB, ICC/IF

**推荐稀释比:** WB: 1/500-1/1000 IF: 1/50-1/200

**种属反应性:** Rabbit

**克隆性:** Rabbit Monoclonal

**克隆编号:** R05-8F8

**分子量:** Calculated MW: 50 kDa; Observed MW: 50 kDa

**亚型:** IgG


**纯化:** Affinity Purified

**种属反应性:** Human, Mouse and Rat


**成分:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**研究领域:** Epigenetics and Nuclear Signaling

**储存和运输:** Store at -20°C. Avoid repeated freezing and thawing



Immunocytochemistry analysis of AdSS 2 (green) in HeLa using AdSS 2 antibody, and DAPI (blue)



Western blot analysis of AdSS 2 in K562, rat Brain, C6, 3T3, HeLa lysates using AdSS 2 antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---