

AF10 RABBIT MAB

货号: N261807

产品全名: AF10 兔单克隆抗体

基因符号 MLLT10; AF10; Protein AF-10; ALL1-fused gene from chromosome 10 蛋白

UNIPROT ID: P55197

背景: Probably involved in transcriptional regulation. In vitro or as fusion protein with KMT2A/MLL1 has transactivation activity. Binds to cruciform DNA. In cells, binding to unmodified histone H3 regulates DOT1L functions including histone H3 'Lys-79' dimethylation (H3K79me2) and gene activation (PubMed:26439302).

抗原: Recombinant protein of human AF10

经过测试的应用: WB,IHC-P

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R07-2G2

分子量: Calculated MW: 113 kDa; Observed MW: 140 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics and Nuclear Signaling

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of AF10 in HeLa, K562 lysates using AF10 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
