

AHA1 RABBIT MAB

货号: N261809

产品全名: AHA1 兔单克隆抗体

基因符号 AHA1; C14orf3; HSPC322; Activator of 90 kDa heat shock protein ATPase homolog 1; AHA1; p38

UNIPROT ID: O95433

背景: Cochaperone that stimulates HSP90 ATPase activity (By similarity). May affect a step in the endoplasmic reticulum to Golgi trafficking.

抗原: A synthetic peptide of human AHA1

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R01-8D9

分子量: Calculated MW: 38 kDa; Observed MW: 42 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of AHA1 in HeLa, HepG2 lysates using AHA1 antibody.



Western blot analysis of AHA1 in rat brain, mouse skin lysates using AHA1 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
