

## ALIX RABBIT PAB

货号: N225934

产品全名: ALIX 兔多抗

基因符号 AIP1; ALIX; HP95; DRIP4

**UNIPROT ID:** Q8WUM4

**背景:** Class E VPS protein involved in concentration and sorting of cargo proteins of the multivesicular body (MVB) for incorporation into intraluminal vesicles (ILVs) that are generated by invagination and scission from the limiting membrane of the endosome.

**抗原:** A synthesized peptide derived from human ALIX

**经过测试的应用:** WB,FC

**推荐稀释比:** WB: 1/500-1/1000 FC: 1/50-1/100

**种属反应性:** Rabbit

**克隆性:** Rabbit Polyclonal

**分子量:** Calculated MW: 96 kDa; Observed MW: 96 kDa

**亚型:** IgG

**纯化:** Affinity Chromatography

**种属反应性:** Human, Mouse and Rat

**成分:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**研究领域:** Cell Biology

**储存和运输:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of ALIX in K562 lysates using ALIX antibody.