

ANGIOTENSINOGEN RABBIT MAB

货号: N261843

产品全名: Angiotensinogen 兔单克隆抗体

基因符号 Angiotensinogen; Serpin A8; AGT; SERPINA8; Des-Asp[1]-angiotensin II; Serine (or cysteine) proteinase inhibitor

UNIPROT ID: P11859

背景: Essential component of the renin-angiotensin system (RAS), a potent regulator of blood pressure, body fluid and electrolyte homeostasis.

抗原: Recombinant protein of mouse Angiotensinogen

经过测试的应用: WB,IHC-P,IP

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R04-8F8

分子量: Calculated MW: 52 kDa; Observed MW: 52 kDa

亚型: IgG

纯化: Affinity Purified

种属反应性: Mouse

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cardiovascular

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Angiotensinogen in mouse liver lysates using Angiotensinogen antibody.