

ANKRD20A1 RABBIT PAB

货号: N225201

产品全名: ANKRD20A1 兔多抗

基因符号 ANKRD20A1; ANKRD20A; Ankyrin repeat domain-containing protein 20A1

UNIPROT ID: Q5TYW2

背景: Located in plasma membrane.

抗原: Synthesized peptide derived from the Internal region of human A20.

经过测试的应用: WB,ELISA

推荐稀释比: WB: 1/500-1/1000 ELISA: 1/10000

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 94 kDa; Observed MW: 94 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of ANKRD20A1 in 3T3 lysates using ANKRD20A1 antibody.