

ANNEXIN V RABBIT PAB

货号: N225793

产品全名: Annexin V 兔多抗

基因符号 PP4; ANX5; ENX2; RPRGL3.

UNIPROT ID: P08758

背景: This protein is an anticoagulant protein that acts as an indirect inhibitor of the thromboplastin-specific complex, which is involved in the blood coagulation cascade.

抗原: A synthesized peptide derived from human Annexin V

经过测试的应用: WB, ICC/IF, FC

推荐稀释比: WB: 1/500-1/1000 IF: 1/50-1/200 FC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 36 kDa; Observed MW: 36 kDa

亚型: IgG


纯化: Affinity Chromatography

种属反应性: Human, Mouse and Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Annexin V in HeLa lysates using Annexin V antibody.