

ANNEXIN VII RABBIT MAB

货号: N261847

产品全名: Annexin VII 兔单克隆抗体

基因符号 SNX; ANX7; SYNEXIN

UNIPROT ID: P20073

背景: Calcium/phospholipid-binding protein which promotes membrane fusion and is involved in exocytosis.

抗原: Recombinant protein of human Annexin VII

经过测试的应用: WB,IHC-P

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R07-1D2

分子量: Calculated MW: 53 kDa; Observed MW: 53 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human,Rat


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Annexin VII antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of Annexin VII in K562, rat Brain lysates using Annexin VII antibody.