

AP2 ALPHA RABBIT MAB

货号: N261850

产品全名: AP2 alpha 兔单克隆抗体

基因符号 AP-2; BOFS; AP2TF; TFAP2; AP-2alpha

UNIPROT ID: P05549

背景: Sequence-specific DNA-binding protein that interacts with inducible viral and cellular enhancer elements to regulate transcription of selected genes. AP-2 factors bind to the consensus sequence 5'-GCCNNNGGC-3' and activate genes involved in a large spectrum of important biological functions including proper eye, face, body wall, limb and neural tube development.

抗原: A synthetic peptide of human AP2 alpha

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R07-7E5

分子量: Calculated MW: 48 kDa; Observed MW: 48 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics and Nuclear Signaling

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of AP2alpha in C6 lysates using AP2 alpha antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
