

## ARF1 RABBIT PAB

货号: N225711

产品全名: ARF1 兔多抗

基因符号 ADP-ribosylation factor 1

**UNIPROT ID:** P84077

**背景:** Stimulate the ADP-ribosyltransferase activity of cholera toxin and play a role in vesicular trafficking as activators of phospholipase D. GTP-binding protein that functions as an allosteric activator of the cholera toxin catalytic subunit, an ADP-ribosyltransferase. Involved in protein trafficking among different compartments.

**抗原:** A synthesized peptide derived from human ARF1

**经过测试的应用:** WB,IHC-P,FC

**推荐稀释比:** WB: 1/500-1/1000 IHC: 1/50-1/100 FC: 1/50-1/100

**种属反应性:** Rabbit

**克隆性:** Rabbit Polyclonal

**分子量:** Calculated MW: 21 kDa; Observed MW: 21 kDa

**亚型:** IgG


**纯化:** Affinity Chromatography

**种属反应性:** Human, Mouse and Rat

**成分:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**研究领域:** Signal Transduction

**储存和运输:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of ARF1 in HEK293 lysates using ARF1 antibody.