

ARFRP1 RABBIT MAB

货号: N261581

产品全名: ARFRP1 兔单克隆抗体

基因符号 ARP; Arp1; ARL18

UNIPROT ID: Q13795

背景: Trans-Golgi-associated GTPase that regulates protein sorting.

抗原: A synthetic peptide of human ARFRP1

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R08-5C5

分子量: Calculated MW: 23 kDa; Observed MW: 23 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of ARFRP1/ARP in C6 lysates using ARFRP1 antibody.