

ARHGAP25 RABBIT MAB

货号: N261869

产品全名: ARHGAP25 兔单克隆抗体

基因符号 Rho GTPase activating protein 25; KAIA0053

UNIPROT ID: P42331

背景: GTPase activator for the Rho-type GTPases by converting them to an inactive GDP-bound state.

抗原: Recombinant protein of human ARHGAP25

经过测试的应用: WB,IHC-P

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R07-6G6

分子量: Calculated MW: 73 kDa; Observed MW: 73 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human,Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of ARHGAP25 in rat brain, Raoms, Jurkat lysates using ARHGAP25 antibody.