

## ARYL HYDROCARBON RECEPTOR RABBIT PAB

货号: N225787

产品全名: Aryl Hydrocarbon Receptor 兔多抗

基因符号 Ah receptor; AhR; Class E basic helix-loop-helix protein 76; bHLHe76; AHR

**UNIPROT ID:** P35869

**背景:** The aryl hydrocarbon receptor (AhR) is a ligand activated transcription factor involved in xenobiotic metabolism, cell cycle regulation, and development in response to both endogenous and environmental signals. Involved in cell-cycle regulation. Likely to play an important role in the development and maturation of many tissues.

**抗原:** A synthesized peptide derived from human AhR

**经过测试的应用:** WB, ICC/IF

**推荐稀释比:** WB: 1/500-1/1000 IF: 1/50-1/200

**种属反应性:** Rabbit

**克隆性:** Rabbit Polyclonal

**分子量:** Calculated MW: 96 kDa; Observed MW: 105 kDa

**亚型:** IgG


**纯化:** Affinity Chromatography

**种属反应性:** Human

**成分:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**研究领域:** Epigenetics and Nuclear Signaling

**储存和运输:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of AhR in HEK293 lysates using Aryl Hydrocarbon Receptor antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---