

ATAXIN 1 RABBIT MAB

货号: N261881

产品全名: Ataxin 1 兔单克隆抗体

基因符号 ATX1; SCA1; D6S504E

UNIPROT ID: P54253

背景: Chromatin-binding factor that repress Notch signaling in the absence of Notch intracellular domain by acting as a CBF1 corepressor. Binds to the HEY promoter and might assist, along with NCOR2, RBPJ-mediated repression. Binds RNA in vitro. May be involved in RNA metabolism (PubMed:21475249). In concert with CIC and ATXN1L, involved in brain development .

抗原: Recombinant protein of human Ataxin 1

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R01-8H1

分子量: Calculated MW: 87 kDa; Observed MW: 105 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Mouse

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics and Nuclear Signaling

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Ataxin 1 in 3T3 lysates using Ataxin 1 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
