

ATG13 RABBIT MAB

货号: N261586

产品全名: ATG13 兔单克隆抗体

基因符号 Autophagy-related protein 13

UNIPROT ID: O75143

背景: Autophagy factor required for autophagosome formation and mitophagy.

抗原: Recombinant protein of human KIAA0652

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R08-8E8

分子量: Calculated MW: 57 kDa; Observed MW: 72 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of ATG13 in K562, C6, HeLa lysates using ATG13 antibody.



Western blot analysis of ATG13 in Raw264.7 lysates using ATG13 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
