

## AVEN RABBIT MAB

货号: N261899

产品全名: AVEN 兔单克隆抗体

基因符号 PDCD12

**UNIPROT ID:** Q9NQS1

背景: Protects against apoptosis mediated by Apaf-1.

抗原: A synthetic peptide of human AVEN

经过测试的应用: WB,IHC-P

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R06-7A9

分子量: Calculated MW: 39 kDa; Observed MW: 50 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat


成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of AVEN in rat Brain, MCF-7, 3T3, HeLa lysates using AVEN antibody.



Immunohistochemistry analysis of paraffin-embedded mouse testis using AVEN antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---