

## AXIN 2 RABBIT PAB

货号: N225809

产品全名: Axin 2 兔多抗

基因符号 Axil; Axin like protein; Axin2; Axis inhibition protein 2; Conductin

**UNIPROT ID:** Q9Y2T1

背景: Axin1 (Axis inhibition protein 1) and Axin2 are multidomain scaffold proteins that negatively regulate Wnt signaling. Axin1 interacts with APC, GSK-3 $\beta$ , Dvl, and  $\beta$ -catenin and promotes the GSK-3 $\beta$ -mediated phosphorylation and subsequent degradation of  $\beta$ -catenin. Down-regulates beta-catenin.

抗原: A synthesized peptide derived from human Axin2

经过测试的应用: WB,IHC-P,ICC/IF

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 94 kDa; Observed MW: 94 kDa

亚型: IgG


纯化: Affinity Chromatography

种属反应性: Human, Mouse and Rat

成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Axin2 in (1) SW480 lysates; (2) PC-12 lysates using Axin 2 antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---