

BECLIN 1 RABBIT MAB

货号: N261689

产品全名: Beclin 1 兔单克隆抗体

基因符号 BECN1

UNIPROT ID: Q14457

背景: Plays a central role in autophagy. Required for the abscission step in cytokinesis. May play a role in antiviral host defense. Protects against infection by a neurovirulent strain of Sindbis virus.

抗原: Recombinant protein of human Beclin 1

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R02-8D5

分子量: Calculated MW: 52 kDa; Observed MW: 60 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Autophagy&Autophagosome

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Beclin 1 in K562, rat Brain, C6, 3T3, Hela lysates using Beclin 1 antibody.