

BNIP3L RABBIT MAB

货号: N261930

产品全名: BNIP3L 兔单克隆抗体

基因符号 Adenovirus E1B19k binding protein B5; BNIP3 α ; BNIP3H; BNIP3L; BNIP3L protein; NIP3 like protein X; NIP3L; NIX

UNIPROT ID: O60238

背景: Induces apoptosis. Interacts with viral and cellular anti-apoptosis proteins. Can overcome the suppressors BCL-2 and BCL-XL, although high levels of BCL-XL expression will inhibit apoptosis. May function as a tumor suppressor. Inhibits apoptosis induced by BNIP3.

抗原: A synthetic peptide of human BNIP3L

经过测试的应用: WB,IHC-P

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R01-0F5

分子量: Calculated MW: 24 kDa; Observed MW: 38 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using BNIP3L antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of BNIP3L in Jurkat lysates using BNIP3L antibody.