

BTK RABBIT MAB

货号: N261939

产品全名: BTK 兔单克隆抗体

基因符号 BTK; AGMX1; ATK; BPK; Tyrosine-protein kinase BTK;
Agammaglobulinaemia tyrosine kinase; ATK; B-cell progenitor kinase; BPK;
Bruton tyrosine kinase

UNIPROT ID: Q06187

背景: Defects in the Bruton tyrosine kinase (BTK) gene cause Agammaglobulinemia. Agammaglobulinemia is an X-linked immunodeficiency characterized by failure to produce mature B lymphocyte cells and associated with a failure of Ig heavy chain rearrangement.

抗原: Recombinant protein of human BTK

经过测试的应用: WB,IHC-P,IP

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R09-8A9

分子量: Calculated MW: 76 kDa; Observed MW: 76 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human tonsil using BTK antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of BTK in K562 lysates using BTK antibody.