

## CALBINDIN RABBIT MAB

货号: N261946

产品全名: Calbindin 兔单克隆抗体

基因符号 Calbindin

**UNIPROT ID:** P12658

背景: Buffers cytosolic calcium. May stimulate a membrane  $\text{Ca}^{2+}$ -ATPase and a 3',5'-cyclic nucleotide phosphodiesterase.

抗原: Recombinant protein of mouse Calbindin

经过测试的应用: WB,IHC-P,IP

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R06-4D6

分子量: Calculated MW: 30 kDa; Observed MW: 30 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat


成分: PBS (without  $\text{Mg}^{2+}$  and  $\text{Ca}^{2+}$ ), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Neuroscience

储存和运输: Store at  $-20^{\circ}\text{C}$ . Avoid repeated freezing and thawing



Western blot analysis of Calbindin in rat Brain lysates using Calbindin antibody.



Immunohistochemistry analysis of paraffin-embedded mouse kidney using Calbindin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.