

CALBINDIN RABBIT PAB

货号: N225792

产品全名: Calbindin 兔多抗

基因符号 avian-type; CAB27; CALB1; Calbindin 1 28kDa; Calbindin; D28K; RTVL H 蛋白

UNIPROT ID: P05937

背景: Buffers cytosolic calcium. May stimulate a membrane Ca(2+)-ATPase and a 3',5'-cyclic nucleotide phosphodiesterase.

抗原: A synthesized peptide derived from human Calbindin

经过测试的应用: WB,IHC-P,ICC/IF

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 30 kDa; Observed MW: 30 kDa

亚型: IgG


纯化: Affinity Chromatography

种属反应性: Human, Mouse and Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Neuroscience

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Calbindin in Human fetal brain lysates using Calbindin antibody.