

CALRETININ RABBIT PAB

货号: N225840

产品全名: Calretinin 兔多抗

基因符号 Calretinin; CR; 29 kDa calbindin

UNIPROT ID: P22676

背景: Calretinin is a calcium-binding protein which is abundant in auditory neurons.

抗原: A synthesized peptide derived from human Calretinin

经过测试的应用: WB,IHC-P,ICC/IF,FC

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 FC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 32 kDa; Observed MW: 29 kDa

亚型: IgG


纯化: Affinity Chromatography

种属反应性: Human, Mouse and Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Tags & Cell Markers

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Calretinin in SH-SY5Y lysates using Calretinin antibody.