

CAV GAMMA 4 RABBIT PAB

货号: N225535

产品全名: Cav gamma 4 兔多抗

基因符号 acng4 antibody; TARP gamma-4 antibody; Transmembrane AMPAR regulatory protein gamma-4 antibody

UNIPROT ID: Q9UBN1

背景: The protein encoded by this gene is a type I transmembrane AMPA receptor regulatory protein (TARP). TARPs regulate both trafficking and channel gating of the AMPA receptors. This gene is part of a functionally diverse eight-member protein subfamily of the PMP-22/EMP/MP20 family and is located in a cluster with two family members, a type II TARP and a calcium channel gamma subunit.

抗原: Synthetic Peptide of L-type Ca⁺⁺ CP γ4

经过测试的应用: WB,IHC-P

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 37 kDa; Observed MW: 37 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Cav gamma 4 in Human Brain, mouse Brain, rat Brain lysates using Cav gamma 4 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
