

## CBX4 RABBIT MAB

货号: N263254

产品全名: CBX4 兔单克隆抗体

基因符号 CBX4; hPc2; NBPI6; PC2

**UNIPROT ID:** O00257

**背景:** E3 SUMO-protein ligase which facilitates SUMO1 conjugation by UBE2I. Component of the Polycomb group (PcG) multiprotein PRC1 complex, a complex required to maintain the transcriptionally repressive state of many genes, including Hox genes, throughout development.

**抗原:** A synthetic peptide of human CBX4

**经过测试的应用:** WB,IP

**推荐稀释比:** WB: 1/500-1/1000 IP: 1/20

**种属反应性:** Rabbit

**克隆性:** Rabbit Monoclonal

**克隆编号:** R08-413

**分子量:** Calculated MW: 61 kDa; Observed MW: 61 kDa

**亚型:** IgG

**纯化:** Affinity Purified

**种属反应性:** Human, Mouse and Rat

**成分:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**研究领域:** Epigenetics and Nuclear Signaling

**储存和运输:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of CBX4 in HeLa, A549, HL-60, U2OS, C6 lysates using CBX4 antibody.

Western blot analysis of CBX4 in 3T3 lysates using CBX4 antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---