

## CD134/OX40L RECEPTER RABBIT PAB

货号: N225945

产品全名: CD134/OX40L Receptor 兔多抗

基因符号 TNFRSF4; TXGPIL; Tumor necrosis factor receptor superfamily member 4; ACT35 antigen; OX40L receptor; TAX transcriptionally-activated glycoprotein 1 receptor; CD antigen CD134

**UNIPROT ID:** P43489

背景: Receptor for TNFSF4/OX40L/GP34.

抗原: A synthesized peptide derived from human CD134

经过测试的应用: WB,FC

推荐稀释比: WB: 1/500-1/1000 FC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 29 kDa; Observed MW: 45 kDa

亚型: IgG


纯化: Affinity Chromatography

种属反应性: Human

成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Immunology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of CD134 in Human fetal spleen lysates using CD134/OX40L Receptor antibody.