

CD168 RABBIT MAB

货号: N261984

产品全名: CD168 兔单克隆抗体

基因符号 HMMR; IHABP; RHAMM; Hyaluronan mediated motility receptor; Intracellular hyaluronic acid-binding protein; Receptor for hyaluronan-mediated motility; CD antigen CD168

UNIPROT ID: O75330

背景: Involved in cell motility. When hyaluronan binds to HMMR, the phosphorylation of a number of proteins, including PTK2/FAK1 occurs. May also be involved in cellular transformation and metastasis formation, and in regulating extracellular-regulated kinase (ERK) activity.

抗原: A synthetic peptide of human CD168

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R02-3E1

分子量: Calculated MW: 84 kDa; Observed MW: 84 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Stem Cells

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of CD168 in K562 lysates using CD168 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
