

CD276 RABBIT MAB

货号: N261991

产品全名: CD276 兔单克隆抗体

基因符号 B7h3; B7RP-2; AU016588; 6030411F23Rik

UNIPROT ID: Q8VE98

背景: Modulates T-cell-mediated immune responses and the development of acute and chronic transplant rejection. Plays a positive regulatory role in bone formation and has a dual role in the bone-immune interface. Induces antitumor immunity as it activates both acquired and innate immunity leading to natural killer cell and CD8 T-cell dependent killing of tumor cells.

抗原: Recombinant protein of mouse CD276

经过测试的应用: WB,IHC-P,IP

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R06-9A3

分子量: Calculated MW: 34 kDa; Observed MW: 90-110 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Immunology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of CD276 in HepG2 lysates using CD276 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
