

## CD3G RABBIT MAB

货号: N262001

产品全名: CD3G 兔单克隆抗体

基因符号 CD3 gamma; CD3g; T3G

**UNIPROT ID:** P09693

背景: The CD3 complex mediates signal transduction.

抗原: A synthetic peptide of human CD3G

经过测试的应用: WB,ICC/IF

推荐稀释比: WB: 1/500-1/1000 IF: 1/50-1/200

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R02-4D1

分子量: Calculated MW: 20 kDa; Observed MW: 20 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human


成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Immunology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of CD3G in Jurkat lysates using CD3G antibody.



Immunocytochemistry analysis of CD3G (green) in Jurkat using CD3G antibody, and DAPI (blue).