

CD42D RABBIT PAB

货号: N225117

产品全名: CD42d 兔多抗

基因符号 GP5; Platelet glycoprotein V; GPV; Glycoprotein 5; CD42d

UNIPROT ID: P40197

背景: The GPIb-V-IX complex functions as the vWF receptor and mediates vWF-dependent platelet adhesion to blood vessels.

抗原: The antiserum was produced against synthesized peptide derived from the Internal region of human GP5. AA range:331-380

经过测试的应用: WB,IHC-P,ELISA

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 ELISA: 1/10000

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 61 kDa; Observed MW: 61 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cardiovascular

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human tonsils using CD42d antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of CD42d in HeLa lysates using CD42d antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
