

## CD73 RABBIT PAB

货号: N225108

产品全名: CD73 兔多抗

基因符号 NT5E; NT5; NTE; 5'-nucleotidase; 5'-NT; Ecto-5'-nucleotidase; CD73

**UNIPROT ID:** P21589

背景: Hydrolyzes extracellular nucleotides into membrane permeable nucleosides.

抗原: The antiserum was produced against synthesized peptide derived from the Internal region of human NT5E. AA range:431-480

经过测试的应用: WB,IHC-P,ELISA

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 ELISA: 1/10000

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 63 kDa; Observed MW: 68 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human


成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cardiovascular

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human skin using CD73 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of CD73 in MCF-7, L929, K562 lysates using CD73 antibody.