

## CD81 RABBIT MAB

货号: N262007

产品全名: CD81 兔单克隆抗体

基因符号 CD81; TAPA1; TSPAN28; CD81 antigen; 26 kDa cell surface protein TAPA-1; Target of the antiproliferative antibody 1; Tetraspanin-28; Tspan-28; CD81

**UNIPROT ID:** P60033

背景: The proteins mediate signal transduction events that play a role in the regulation of cell development, activation, growth and motility. Many research studies demonstrate a role for CD81 in lymphocyte signaling. CD81 is also a well-characterized receptor for Hepatitis C Virus and facilitates the entry of the virus into target cells.

抗原: Recombinant protein of human CD81

经过测试的应用: WB, IHC-F, IHC-P, ICC/IF

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R02-3D6

分子量: Calculated MW: 26 kDa; Observed MW: 26 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse

成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Tags & Cell Markers

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of CD81 in C2C12, 3T3, MCF-7 lysates using CD81 antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---