

CD96 RABBIT PAB

货号: N225119

产品全名: CD96 兔多抗

基因符号 CD96; T-cell surface protein tactile; Cell surface antigen CD96; T cell-activated increased late expression protein; CD96

UNIPROT ID: P40200

背景: May be involved in adhesive interactions of activated T and NK cells during the late phase of the immune response.

抗原: The antiserum was produced against synthesized peptide derived from the Internal region of human CD96. AA range:291-340

经过测试的应用: WB,ELISA

推荐稀释比: WB: 1/500-1/1000 ELISA: 1/10000

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 66 kDa; Observed MW: 66 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Immunology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of CD96 in 293UV lysates using CD96 antibody.



Western blot analysis of CD96 in 293UV lysates using CD96 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
