

## CDK4 RABBIT MAB

货号: N262018

产品全名: CDK4 兔单克隆抗体

基因符号 CMM3; PSK-J3

**UNIPROT ID:** P11802

**背景:** Cyclin-dependent kinase activity is regulated by T-loop phosphorylation (Thr172 in the case of CDK4), by the abundance of their cyclin partners, and by association with CDK inhibitors of the Cip/Kip or INK family of proteins. Cyclin D-CDK4 complexes are major integrators of various mitogenic and antimitogenic signals.

**抗原:** A synthetic peptide of human Cdk4

**经过测试的应用:** WB, ICC/IF

**推荐稀释比:** WB: 1/500-1/1000 IF: 1/50-1/200

**种属反应性:** Rabbit

**克隆性:** Rabbit Monoclonal

**克隆编号:** R07-4G7

**分子量:** Calculated MW: 34 kDa; Observed MW: 34 kDa

**亚型:** IgG


**纯化:** Affinity Purified

**种属反应性:** Human

**成分:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**研究领域:** Cell Biology

**储存和运输:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Cdk4 in HeLa lysates using CDK4 antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---