

CDKN2A/P16INK4A RABBIT PAB

货号: N226007

产品全名: CDKN2A/p16INK4a 兔多抗

基因符号 CDKN2A; CDKN2; MTS1; Cyclin-dependent kinase inhibitor 2A; isoforms 1/2/3; Cyclin-dependent kinase 4 inhibitor A; CDK4I; Multiple tumor suppressor 1; MTS-1; p16-INK4a; p16-INK4; p16INK4A

UNIPROT ID: P42771

背景: p16-INK4A a cell-cycle regulatory protein that Interacts with CDK4 and CDK6, inhibiting their ability to interact with cyclins D. Inhibits the phosphorylation of the retinoblastoma protein by CDK4 or CDK6. Four alternatively spliced isoforms have been reported.

抗原: A synthetic peptide of human CDKN2A/p16INK4a

经过测试的应用: FC,ICC/IF,IP,WB

推荐稀释比: WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20 FC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 17 kDa; Observed MW: 17 kDa

亚型: IgG

纯化: Affinity Purified

种属反应性: Human


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of p16 INK4A in 293 lysates using p16 INK4A antibody.



Immunocytochemistry analysis of CDKN2A (green) in 293 using CDKN2A antibody, and DAPI (blue).



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
