

CHROMOGRANIN A RABBIT MAB

货号: N262036

产品全名: Chromogranin A 兔单克隆抗体

基因符号 beta Granin; Chromogranin A; ChrA; SP-I; Vasostatin I;
Vasostatin-2; Vasostatin II; Pancreastatin; CHGA

UNIPROT ID: P10645

背景: Pancreastatin strongly inhibits glucose induced insulin release from the pancreas.

抗原: A synthetic peptide of human Chromogranin A

经过测试的应用: WB,IHC-F,IHC-P,ICC/IF

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R01-7E1

分子量: Calculated MW: 51 kDa; Observed MW: 92 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Neuroscience

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Chromogranin A in PC-12 lysates using Chromogranin A antibody.