

COLLAGEN VI ALPHA 1 RABBIT MAB

货号: N262053

产品全名: Collagen VI alpha 1 兔单克隆抗体

基因符号 OPLL; BTHLM1; UCHMD1

UNIPROT ID: P12109

背景: Collagen VI acts as a cell-binding protein.

抗原: Recombinant protein of human Collagen VI

经过测试的应用: WB,IHC-F,IHC-P,ICC/IF

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R02-819

分子量: Calculated MW: 109 kDa; Observed MW: 147 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Collagen VI in Human placenta lysates using Collagen VI alpha 1 antibody.