

## COLLAGEN XIX ALPHA 1 RABBIT PAB

货号: N225231

产品全名: Collagen XIX alpha 1 兔多抗

基因符号 Collagen alpha-1(XIX) chain (Collagen alpha-1(Y) chain)

**UNIPROT ID:** Q14993

背景: May act as a cross-bridge between fibrils and other extracellular matrix molecules.

抗原: Synthesized peptide derived from human COL19A1. at AA range: 161-210

经过测试的应用: WB,ELISA

推荐稀释比: WB: 1/500-1/1000 ELISA: 1/10000

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 115 kDa; Observed MW: 130 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Collagen XIX alpha 1 in HELA lysates using Collagen XIX alpha 1 antibody.