

COLLAGEN XXVIII ALPHA 1 RABBIT PAB

货号: N225209

产品全名: Collagen XXVIII alpha 1 兔多抗

基因符号 Collagen alpha-1(XXVIII) chain

UNIPROT ID: Q2UY09

背景: May act as a cell-binding protein.

抗原: The antiserum was produced against synthesized peptide derived from the N-terminal region of human COL28A1. AA range:111-160

经过测试的应用: WB,ELISA

推荐稀释比: WB: 1/500-1/1000 ELISA: 1/10000

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 117 kDa; Observed MW: 117 kDa

亚型: IgG

纯化: Affinity Purified

种属反应性: Human,Mouse

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Collagen XXVIII alpha 1 in 3T3 HEPG2 lysates using Collagen XXVIII alpha 1 antibody.